



Designation: E1064 – 16

Standard Test Method for Water in Organic Liquids by Coulometric Karl Fischer Titration¹

This standard is issued under the fixed designation E1064; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope*

1.1 This test method covers the determination of water from 0 to 2.0 % mass in most liquid organic chemicals, with Karl Fischer reagent, using an automated coulometric titration procedure. Use of this test method is not applicable for liquefied gas products such as Liquid Petroleum Gas (LPG), Butane, Propane, Liquid Natural Gas (LNG), etc.

1.2 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

1.3 Review the current Safety Data Sheets (SDS) for detailed information concerning toxicity, first-aid procedures, handling, and safety precautions.

1.4 *This standard does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety, health, and environmental practices and determine the applicability of regulatory limitations prior to use.* Specific precautionary statements are given in Section 8.

1.5 *This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.*

2. Referenced Documents

2.1 *ASTM Standards:*²

D1193 Specification for Reagent Water

D4672 Test Method for Polyurethane Raw Materials: Determination of Water Content of Polyols

¹ This test method is under the jurisdiction of ASTM Committee D16 on Aromatic, Industrial, Specialty and Related Chemicals and is the direct responsibility of Subcommittee D16.04 on Instrumental Analysis.

Current edition approved April 1, 2016. Published May 2016. Originally approved in 1985. Last previous edition approved in 2012 as E1064 – 12. DOI: 10.1520/E1064-16.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

E180 Practice for Determining the Precision of ASTM Methods for Analysis and Testing of Industrial and Specialty Chemicals (Withdrawn 2009)³

E203 Test Method for Water Using Volumetric Karl Fischer Titration

3. Summary of Test Method

3.1 This test method is based on the Karl Fischer reaction for determining water—the reduction of iodine by sulfur dioxide in the presence of water to form sulfur trioxide and hydriodic acid. The reaction becomes quantitative only when pyridine or other organic base and methanol or other alcohol are present. Unlike the volumetric Karl Fischer reagents that include iodine, the coulometric technique electrolytically generates iodine, with 10.71 C of generating current corresponding to 1 mg of water in accordance with Faraday's law.

4. Significance and Use

4.1 The coulometric technique is especially suited for determining low concentrations of water in organic liquids that would yield small titers by the Karl Fischer volumetric procedure. The precision and accuracy of the coulometric technique decreases for concentrations of water much greater than 2.0 % because of the difficulty in measuring the small size of sample required. The test method assumes 100 % efficiency of coulombs in iodine production. Provision is made for verifying this efficiency. (See Table 1 and Note 5.)

5. Interferences

5.1 Interfering substances are the same as are encountered in the volumetric Karl Fischer titration. A detailed discussion of interfering substances can be found in the treatise on "Aquametry."⁴

5.2 Test Method E203 discusses organic compounds in which water may be determined directly and compounds in

³ The last approved version of this historical standard is referenced on www.astm.org.

⁴ J. Mitchell, Jr. and D. M. Smith, "Aquametry"—*A Treatise on Methods for the Determination of Water, Part III*—The Karl Fischer Reagent, 2nd Ed., J. Wiley and Sons, Inc., New York, NY 1990.

*A Summary of Changes section appears at the end of this standard

which water cannot be determined directly, but in which interferences may be eliminated by suitable chemical reactions.

6. Apparatus

6.1 *Automatic Titrator*,⁵ consisting of a control unit, titration vessel, dual platinum sensing electrode, generator assembly, and magnetic stirrer. The instrument is designed to coulometrically generate iodine that reacts stoichiometrically with the water present in the sample solution. The coulombs of electricity required to generate the reagent are converted to micrograms of water, which is obtained as a direct digital readout.

6.2 *Syringe*, 50-mL, fitted with an 115-mm hypodermic needle for removing excess solution from the titration chamber.

NOTE 1—Rinse all glass syringes and needles with anhydrous acetone after cleaning, then dry in an oven at 100°C for at least 1 h and store in a desiccator. Plastic syringes shall be disposed of following use.

6.3 *Syringe*, 20-mL, fitted with an 115-mm hypodermic needle for introduction of neutralizing solution into the titration chamber (see [Note 1](#)).

6.4 *Syringes*, 1- and 5-mL, fitted with 115-mm hypodermic needles for introduction of samples into titration chamber (see [Note 1](#)).

6.5 *Syringe*, 5 μ L, fitted with 115-mm hypodermic needle for standardization of instrument (see [Note 1](#)).

6.6 *Fluorocarbon Sealing Grease or TFE-Fluorocarbon*, to seal the titration chamber against atmospheric moisture.

6.7 *Septa*, to seal sample port but allow introduction of samples by a needle with a minimum of moisture contamination. Replace serum caps and septa as required to prevent air leakage as indicated by instrument drift.

6.8 *Serum Bottles*.

6.9 *Oven*, temperature $100 \pm 5^\circ\text{C}$.

6.10 *Dessicator*, standard laboratory type with color change indicator.

6.11 *Analytical Balance*, capable of weighing to ± 0.0001 g.

7. Reagents

7.1 *Purity of Reagents*—Unless otherwise indicated, it is intended that all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available.⁶ Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

⁵ Basic references to the automatic coulometric titrator: M. T. Kelley, R. W. Stelzner, W. R. Laing, and D. J. Fisher, *Analytical Chemistry* 31, No. 2, 220 (1959) and A. W. Meyer, Jr. and C. M. Boyd, *Analytical Chemistry* 31, No. 2, 215 (1959).

⁶ *Reagent Chemicals, American Chemical Society Specifications*, American Chemical Society, Washington, DC. For suggestions on the testing of reagents not listed by the American Chemical Society, see *Analar Standards for Laboratory Chemicals*, BDH Ltd., Poole, Dorset, U.K., and the *United States Pharmacopeia and National Formulary*, U.S. Pharmacopeial Convention, Inc. (USP), Rockville, MD.

7.2 *Purity of Water*—Unless otherwise indicated, reference to water shall be understood to mean Type II or Type III reagent water, conforming to Specification [D1193](#).

7.3 *Karl Fischer Reagents*—Commercial coulometric KF reagents and reagent systems of various types are available for use with autotitrators for water determination. Traditionally, pyridine was the organic base used in KF reagents. Pyridine-free formulations are available and are preferred by most KF instrument manufacturers for use with their equipment. The pyridine-free reagents are less toxic, less odorous, and more stable than those containing pyridine. The use of pyridine-free reagents is recommended whenever possible. Coulometric titrations normally require two reagent solutions. An anolyte or solvent titration solution and a catholyte or generator titrant solution. However, now reagents can be purchased in one or two component reagent systems. A one component reagent system contains all the components required for a Karl Fischer titration in a single solution. A two component system incorporates separate solutions for the solvent and titrant.

NOTE 2—Two good references on pyridine-free reagents are the Hydranal[®] Manual-Eugen Schotz Reagents for Karl Fischer Titration, from Riedel-deHaen (www.rdhlab.de) or Sigma Aldrich (www.sigma-aldrich.com) and Moisture Measurement by Karl Fischer Titrimetry, 2nd ed., by GFS Chemicals, Inc., January 2004.

7.3.1 *Generator Titrant Solution (catholyte)*, containing iodine, sulfur dioxide, pyridine or other organic base and methanol or other alcohol to provide iodine in the reaction mixture.

7.3.2 *Solvent Titration Solution (anolyte)*, prepared as per instrument specifications.

7.3.3 *Neutralizing Solution*, methanol containing approximately 20 mg $\text{H}_2\text{O}/\text{mL}$.

8. Safety Precautions

8.1 The reagents contain one or more of the following: iodine, pyridine or other organic base, sulfur dioxide, and methanol or other alcohol. Wear chemical resistant gloves when mixing the reagents and removing solution from the titration chamber. Care must be exercised to avoid inhalation of reagent vapors, or direct contact of the reagent with the skin.

9. Sampling

9.1 Because of the low concentration of water to be measured, maximum care must be exercised at all times to avoid contaminating the sample with moisture from the sample container, the atmosphere, or transfer equipment.

9.1.1 Dry the sample bottles and caps overnight in an oven at 100°C before using. Allow to cool in a desiccator before filling and sealing.

9.1.2 Fill the sample bottle as rapidly as possible to within 25 mm of the top and immediately seal.

9.2 Remove the test specimens for analysis from the sample bottle with a dry hypodermic syringe. Inject dry nitrogen into the sample bottle with the syringe to displace the sample that is removed.

10. Preparation of Apparatus

10.1 Clean, dry, and assemble the titration chamber as directed in the manufacturer's instructions. Use fluorocarbon

grease or TFE-fluorocarbon to seal the upper and lower sections of the chamber from atmospheric moisture. Connect the leads from the sensing and generator electrodes to the titrator.

10.2 Prepare the titration solution (7.3.2) as directed by the instrument manufacturer and fill the instrument reservoir as directed by the manufacturer.

10.3 Add the generator solution (7.3.1) to the generator assembly and replace the cover cap. The surface of the generator solution must be below the level of the titration solution to prevent backflow contamination of the titration solution. The generator solution may have to be removed periodically to maintain the lower level.

10.4 Follow the manufacturer's procedure to dry the titration cell.

11. Verification of Calibration

11.1 Different autotitrators may vary in calibration procedures. Consult the operating manual for the autotitrator in use. Stable, prepackaged Quality Control (QC) water standards are commercially available with 10 mg/kg, 100 mg/kg and 1 % (m/m) water content for this purpose. It is desirable to verify calibration with a standard solution that approximates the same range of water expected to be in the samples.

11.2 It is recommended that a control chart measuring a QC standard sample be established and maintained according to generally accepted guidelines.⁷ Measure the control sample each time a test sample(s) is tested. If the measured value exceeds ± 5 % of the known amount, take appropriate action before proceeding with the sample test.

NOTE 3—This may require replacing or regenerating the reagent solutions.

12. Procedure

12.1 Assemble a dry syringe and needle and attach a small piece of rubber septum to the needle tip. Withdraw 1 to 2 mL of the sample into the syringe and discard the contents into a waste container. Using the following table as a guide, withdraw the proper amount of test sample into the syringe and seal with the rubber septum. Obtain a tare weight to ± 0.1 mg. See Table 1, Note 4 and Note 5.

NOTE 4—If multiple analyses are performed on the same test sample or if an appreciable volume of test sample is withdrawn, a slight pressure should be maintained on the sample bottle by means of a hypodermic needle attached to a dry nitrogen source.

NOTE 5—Alternatively, if syringes of suitable accuracy are available, a

measured volume of sample can be injected and the mass calculated from the volume and density.

12.2 With the analyzer stabilized, carefully insert the needle of the sample syringe through the septum and below the level of solution in the titration chamber. Inject the sample slowly into the titration solution and begin titration. Withdraw the syringe needle, seal and again weigh to the nearest ± 0.1 mg to determine the exact sample mass. Allow the titration to proceed until the end-point is indicated. Record the micrograms of water determined.

NOTE 6—After numerous analyses, the level of solvent accumulated in the titration chamber may have to be reduced. This can be accomplished with a 50-mL syringe or by partially draining the solution if the titration chamber is provided with a stopcock. Discard the solution and replace with fresh titration solution if a stable reading cannot be obtained.

NOTE 7—Replace the generator solution when it becomes yellow and a stable reading cannot be obtained.

13. Calculation

13.1 Calculate the water content of the sample to the nearest 0.001 % mass as follows:

$$\text{water content, \% mass} = \frac{\mu\text{g water found}}{\text{g sample} \times 10000} \quad (1)$$

14. Report

14.1 Report the percentage of water to the nearest 0.001 % mass.

15. Precision and Bias

15.1 *Precision*—The following criteria should be used to judge the acceptability of results when using pyridine-based reagents (see Note 8).

15.1.1 *Repeatability (Single Analyst)*—The coefficient of variation for a single determination has been estimated to be 1.39 % relative at 61 DF. The 95 % limit for the difference between two such runs is 3.9 % relative.

15.1.2 *Laboratory Precision (Within-Laboratory, Between Days Variability)*—The coefficient of variation of results (each the average of duplicate determinations), obtained by the same analyst on different days, was estimated to be 2.00 % relative at 31 DF. The 95 % limit for the difference between two such averages is 5.6 % relative.

15.1.3 *Reproducibility (Multilaboratory)*—The coefficient of variation of results (each the average of duplicate determinations), obtained by analysts in different laboratories, has been estimated to be 6.12 % relative at 7 DF. The 95 % limit for the difference between two such averages is 17.1 % relative.

15.2 *Bias*—The bias of this test method has not been determined since no acceptable reference material has been identified.

NOTE 8—The above precision estimates are based on an inter-laboratory study on samples of toluene, ethyl acetate, ethanol, and acetonitrile, containing approximately 0.02, 0.16, 0.52, and 2.04 % mass water, respectively. One analyst in each of 8 laboratories performed duplicate determinations on two different days, for a total of 128 determinations. Five values for toluene were judged to be outliers. Practice E180 and a special computer program for the analysis of variance of unbalanced nested experiment designs were used to develop these precision estimates.

TABLE 1 Sample Size Estimation

Expected Water Content	Sample Size, mL
0 to 500 mg/kg	5
500 to 1000 mg/kg	2
1000 to 2000 mg/kg	1
0.2 to 0.5 % mass	0.5
0.5 to 2.0 % mass	0.1

⁷ ASTM Manual on Presentation of Data and Control Chart Analysis, 7th Edition, ASTM Manual Series MNL 7A, (revision of Special Technical Publication (STP) 15D).

15.3 *Precision*—The following is an example of the precision attained in an interlaboratory study for determining water with pyridine-free reagents (see [Note 9](#)).

15.3.1 *Repeatability*—Duplicate results obtained by the same analyst should be considered suspect if they differ by more than the percent relative listed in [Table 2](#) for the water content which most closely matches the sample be analyzed.

TABLE 2 Repeatability and Reproducibility Statistics

Water Content, % mass	95 % Limit	
	Repeatability, % relative	Reproducibility, % relative
0.03	2.8	15.9
0.42	3.1	3.2
1.6	3.1	5.2

15.3.2 *Reproducibility*—The average result of duplicates obtained in one laboratory should not differ from that of another laboratory by more than the relative percentage given in [Table 2](#) for the water levels listed.

NOTE 9—The above precision estimates are from Test Method B of Test Method [D4672](#). An interlaboratory study was done in 2000 on three polyol samples containing approximately 0.03, 0.42, and 1.6 % mass water. One analyst performed duplicate determinations and repeated them on a second day. Practice [E180](#) was used to develop the precision estimates.

15.3.3 *Bias*—The bias of this test method has not been determined.

16. Keywords

16.1 coulometric; Karl Fischer; organic liquids; pyridine; pyridine-free; water

APPENDIX

(Nonmandatory Information)

X1. PRECISION AND BIAS FOR GLYCOLS

X1.1 In 2007, ASTM International Committee E15 on Industrial and Specialty Chemicals conducted and completed Interlaboratory Study No. 52 to determine precision data for six test methods used in the analysis of glycols. The precision of this test method is based on the interlaboratory study of E1064, Standard Test Method for the Determination of Water in Organic Liquids by Coulometric Karl Fischer Titration, conducted in 2007. Each of seventeen laboratories were asked to test three different materials. Fourteen laboratories tested MEG, 13 laboratories tested DEG and 13 laboratories tested TEG. Every “test result” represents an individual determination. Two test results were conducted on each of two days for a total of four test results per assay. Note that in the combined study, 8 laboratories used a single analyst, 7 laboratories used two analysts (on different days) and 2 laboratories did not record this information. In the event that there were missing values for one or more laboratories, this information was noted in the results.⁸

X1.1.1 *Repeatability*—Two test results obtained within one laboratory shall be judged not equivalent if they differ by more than the “r” value for that material; “r” is the interval representing the critical difference between two test results for

the same material, obtained by the same operator using the same equipment on the same day in the same laboratory.

X1.1.2 *Reproducibility*—Two test results shall be judged not equivalent if they differ by more than the “R” value for that material; “R” is the interval representing the difference between two test results for the same material, obtained by different operators using different equipment in different laboratories.

X1.1.3 *Intermediate Precision*—The day-to-day standard deviation within a laboratory for results produced by the same operator, determined through statistical analysis following ASTM [E180](#), Standard Practice for Determining the Precision of ASTM Methods for Analysis and Testing of Industrial and Specialty Chemicals. Practice [E180](#) was used to conform to this particular study design which required an estimate of intermediate precision. The statistical analysis was conducted using the SAS statistical analysis software, Version 8.0.

X1.1.3.1 The Practice [E180](#) analysis considers the two test results from each day as being run under repeatability, intermediate, and reproducibility precision for each assay. The repeatability precision would be estimated from the two sets of duplicate test results within each day, and the intermediate precision would be estimated from the agreement between the two days, all pooled over laboratories. Caveat: Since two days is a short time period, the intermediate precision would probably be underestimated by the [E180](#) analysis.

⁸ Supporting data have been filed at ASTM International Headquarters and may be obtained by requesting Research Report No. RR:E15-1062. Contact ASTM Customer Service at service@astm.org.

TABLE X1.1 E1064 Water in Organic Liquids by Coulometric Karl Fischer Titration

Test Result % mass	Sample	Average over all Laboratories	Repeatability Standard Deviation	Intermediate Standard Deviation	Reproducibility Standard Deviation	Repeatability Limit	Intermediate Unit	Reproducibility Limit
Water	MEG	0.0086	0.0009	0.0014	0.0025	0.0026	0.0038	0.0071
Water	DEG	0.0649	0.0012	0.0014	0.0049	0.0032	0.0039	0.0137
Water	TEG	0.0498	0.0019	0.0129	0.0157	0.0054	0.0361	0.0439

X1.1.4 Any judgment in accordance with these two statements would have an approximate 95 % probability of being correct.

X1.2 Bias—At the time of the study, there was no accepted reference material suitable for determining the bias for this test method, therefore no statement on bias is being made.

X1.3 The precision statement was determined through statistical examination of qualified results, from seventeen

laboratories, on three materials. These three materials were described as the following:

Fluid 1:	Monoethylene Glycol
Fluid 2:	Diethylene Glycol
Fluid 3:	Triethylene Glycol

X1.3.1 To judge the equivalency of two test results, it is recommended to choose the material closest in characteristics to the test material.

SUMMARY OF CHANGES

Committee E15 has identified the location of selected changes to this standard since the last issue (E1064-12) that may impact the use of this standard. (Approved April 1, 2016.)

(1) Revised 10.4.

(2) Deleted Note 3, 10.5, and 10.6

(3) Moved Table 1 to Section 12 Procedure.

ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, Tel: (978) 646-2600; <http://www.copyright.com/>